Floc-Tester

AL40 • AL50

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GB Instruction Manual
Page 17–29

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ES Instrucciones
Pagina 45–56

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- **General Notes** 18
 - Safety Instructions 19
 - Introduction 19
 - Instrument description 19
- **Functional description** 20
- **Assembly and installation** 21
 - Electric connection 21
 - Start-up 21
- **Operating controls** 22
- **User notes** 23
- **Technical data** 24
- **Analytical methods** 25
 - Coagulation and flocculation of wastewaters 25
 - Laboratory evaluation by jar test 26
 - Leaching test 27
- **Calculation of velocity gradient (G)** 28
- **Declaration of conformity** 29

Safety Instructions

1. Before using the unit, please read the operating manual supplied with the apparatus carefully.
2. Do not dispose of this equipment as urban waste
3. In order to prevent possible risk of electric shocks, fire and personal injury when the unit is being used, basic safety measures must always be taken, including:
4. Do not turn the instrument on before positioning the beakers for protecting the blades. Ensure that the power switch is turned to "0" (OFF) before removing the beakers containing the stirring liquid.
5. Check that the rating of external power supply corresponds to the rating of the electric line.
6. Do not use the unit if it is not working correctly. In case of malfunctioning, contact your nearest assistance centre.
7. Personal protective equipment must be compatible with the possible risks posed by the material being processed and the glass parts.
8. Follow the cleaning instructions described in this manual.
9. This unit must only be used for laboratory applications
10. The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.
11. This unit has been designed and produced in compliance with the following standards.
12. Safety requirements for electrical apparatus for:

Measurement and control and for laboratory use	CEI EN 61010-1
Electrical equipment for laboratory use	UL 3101-1
General requirement – Canadian electrical code	CAN/CSA-C22.2

Note

- The manufacturer is committed to constantly improving the quality of the products and reserves the right to modify the characteristics without prior notice.

Introduction

The equipment is aimed to satisfy two different analytical operations:

Jar tests for the optimization of dosing coagulants and polyelectrolytes in wastewaters treatment plants.

Leaching tests using diluted acetic acid solutions or water saturated by carbon dioxide for the evaluation of leacheability of heavy metals by solid wastes to be sent to different kinds of dumps.

The stirring speed of the ET flocculators is reproducible meaning that standard test conditioning can be used on a regular basis, an indispensable requisite for obtaining reproducible results.

The examined samples can be lighted by a back fluorescent lamp controlled by a switch located on front panel (1).

The instruments have an ergonomic design, the control pannel is slightly inclined for easy use and good visibility.

The programmable speed of rotation from 10 to 300 r.p.m. (within 1 r.p.m.) is continuously compared to the set value and shown by a display. The unit is equipped by an electronic timer allowing to program the stirring time on two different scales (minutes or hours) or to choose a continuous operation.

Instrument description

The flocculator is manufactured in a strong metallic frame with an epoxy painting which has been rightly studied to give a high degree of protection against chemical and mechanical corrosion to the instrument.

The stainless steel stirring rods are adjustable in height by a self blocking chuck.

Rotation speed is controlled by a gear motor powered by direct current controlled electronically by a microprocessor.

From the frontal panel it's possible to select the same stirring speed from 10 to 300 rpm for each position, the real speed is showed on the relative display.

The electronic timer allow to set three different running time: hour, minutes or continuous way.

General Switch

The general switch allows to turn the instrument on and off. If the switch is set to "OFF", the instrument is off; if the switch is set to "ON", the instrument is on.

Pannel Lighting Switch

This switch allows to turn on the back lamp.

Electronic Timer Window

Trough the keys and the displays of this window (Set ↑ e ↓) it's possible to select the following modality of running:

- Continuous
- With time base. Up to 999 minutes
- With hours base. Up to 99 hours

Stirrer Speed Window

It's possible select the stirrer speed for the stirring rods by the keys and display of this window (↑ e ↓). It is possible to select speed from 10 to 300 rpm. The display shows the stirring rods speed.

Information regarding construction materials

Structure	Steel
Frontal panel label	PET
Foot	NBR

If any liquid being processed should fall onto the above materials, clean them immediately to avoid chemical corrosion.

Electric connection

After unpacking the instrument, place the unit on a laboratory bench.

Before connecting the instrument to the electric line, make sure that the values on the rating plate correspond to those of the electric line.

Use the supplied external power to connect the instrument to the socket.

The supplied external power allow to connect the instrument to net power lines between 100 and 240 Volt with 50 or 60 Hz.

Start-up

After connecting the instrument to the electric line, place the sample containing beakers on the instrument, with the stirring rod facing inside.

Do not turn the instrument on before positioning the beakers (thus protecting the blades).

In order to place the beakers it is necessary to pull stirring rods up, just sufficiently, by holding them by the suitable black pommels; to place the beakers (one for each stirring rod) and to pull down the rods.

Ensure that the power switch is turned to "0" (OFF) before removing the beakers containing the stirring liquid.

Select the the desired modality of running (continuous, with minutes base and with hours base) by the keys of the electronic power window on the control pannel.

Select the desired speed by the keys of stirrer speed window.

The start-up is empowered by the keys ↑ and ↓. It's possible to select the speed from 0 to 300 rpm by the speed arrow keys. The display shows the selected speed of the stirring rods. The progression speed is increased by keeping one of the two keys pushed, in order to reach the desired value faster.

Electronic timer setting

It's possible to select three different type of running by the Set key on the "Electronic Timer Window":

Continuous

With timer (minutes base)

With timer (hours base)

It's possible to select the desired modality just by pressing the Set key and the display will show as follows:

Continuous way: "....."

With timer (minutes): "000"

With timer (hours): "h00"

In order to change the timer modality maintain the Set key pushed for a few seconds. To set the desired time use the keys ↑ and ↓.

The display shows the count down and simultaneously the third decimal point turns on to indicate that the count is operative. At the end of time the instrument stops the stirring, the acoustic signal warns the end of the time and the display shows "End".

Speed setting

By the arrows key on the "Speed" window it's possible to select the desired speed from 10 to 300 rpm with 1 rpm resolution. Maintaining pressed one arrow, increase the progression speed in order to reach more quickly the desired value.

The speed is not active if the time is not selected

Black-out

The Floc-Tester has an inside memory that allows to keep the memory of the stirring speed and of the residual time in case of black out, so that the instrument will come to the same sets as soon as the black out ends.

End of work cycle operations

At the end of the work cycle select 0 rpm and turn OFF the general switch.

If the instrument is not used for a long period it's better to disconnect the external power supply from the electric line.

Maintenance

There is no need of routine nor extraordinary maintenance, but periodically cleaning of the unit as described in this manual.

In compliance with the product guarantee law, the reparings of our units must be carried out in our factory, unless previously otherwise agreements with local distributors.

Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable and non-aggressive detergent to clean the unit.

Disposal

The final disposal of the unit or of its components must be carried out in compliance with the directives and laws in force in the country concerned.

GB Technical data

General

Power supply	V	DC 19 ± 0,5
Power	W	19
Dimensions ET 740	mm (bxhxp)	655x404x296
Dimensions ET 750	mm (bxhxp)	935x404x296
Weigh ET 740	Kg	13
Weigh Et 750	Kg	17
Programmable speed	rpm (1/min)	from 10 to 300
Speed selection	rpm (1/min)	1
Programmable time	min / hours	0 ÷ 999 / 0 ÷ 99 or continuous
Time selection	min / hours	1 / 1
Construction material		Stell
Environmental temp. range	°C	+5...+40
Max storage admitted	°C	-10...+60
Max. humidity	%	Max 80
Operating mode		Continuous
Pollution degree CEI EN61010-1		2
Installation category CEI EN61010-1		2

External Power Supply

Inlet	AC 100÷240V ; 50-60 Hz ; 1.5A
Outlet	DC 19V ; 2000 mA

Coagulation and flocculation of wastewaters

Chemical coagulants are added to wastewaters for the primary purpose of removing suspended solids and also phosphorus. The addition of chemicals reduces also heavy metal concentrations and improves disinfection efficiency. Chemical coagulation can be obtained by adding to wastewaters defined amounts of lime (calcium hydroxide), alum (aluminium sulfate) or iron salts (ferric or ferrous).

The coagulation-sedimentation process typically involves:

injection and mixing of the coagulant that neutralize the electric charge, prevalently negative, on suspended particles. Good and rapid mixing of the coagulant and wastewater is important to ensure efficient use of the chemical. Typical detention times in plant basins for turbulent mixing are 15 to 120 seconds.

agglomeration of the coagulated particles into large settable flocs. This agglomeration is accomplished by stirring the water slowly so to allow the solids formed by the addition of chemicals to grow in size so that they can be removed by gravity settling.

The stirring is usually performed by slowly rotating paddles during 10 to 30 minutes.

sedimentation of the flocculated material by gravity separation. The solids are allowed to settle by gravity from the chemically treated standing wastewater. The effluent from the chemical settling tank is very clear and does not usually require further filtration.

The quantity of chemical coagulant required to achieve good coagulation varies with time and

from wastewater to wastewater.

Typical coagulant doses are:

75 to 250 g/m³ for alum,

45 to 90 g/m³ for ferric chloride,

200 to 400 g/m³ for lime.

Laboratory evaluation by Jar test

The choice and dosing of the chemical coagulant to be adopted for the removal of suspended solids from wastewaters derive from the results of laboratory evaluations by the so called Jar test. Multiple stirrers with reproducible stirring speeds allow to adopt standard conditions for the test, that are the basic requirement for reproducible results.

The standard conditions relate to:

- volume of wastewater in each stirring position
- dimensions and shape of containing vessel and stirring paddle
- time and speed of stirring during turbulent mixing phase
- time and speed of stirring during flocculation phase
- evaluation criteria for the obtained results.

The most diffused adopted conditions are:

- 1000 ml glass beakers (jars), tall form, 105 mm of diameter
- 600 ml wastewater samples
- stirring paddles 25 mm high, 75 mm wide, with a thickness of 1 mm
- stirring height of paddles roughly at middle height of wastewater sample
- quick turbulent stirring after chemical addition: 120 rpm for 120 seconds
- slow speed stirring during flocculation: 30 rpm for 25 minutes
- first evaluation of results after 5 minutes of sedimentation after the end of stirring.

The results can be evaluated on the basis of different criteria:

- a) evaluation of floc dimensions with a numerical degree of merit (0 = no flocs, 2 = very small hardly visible, 4 = small flocs, 6 = medium size flocs, 8 = good size flocs, 10 = very large flocs).
- b) time from the addition of chemicals to first appearing of flocs
- c) evaluation of residual turbidity of supernatant, after a determined sedimentation time, by turbidimeter
- d) measurement of electrocyntetic potential of suspended particles on a sample taken immediately after the addition and mixing of chemicals. Coagulations should be easier when zeta potential of particles is very low.
- e) evaluation of filterability of clarified water by standardized membrane filters under pressure. The reduction of water flow with time is related to the degree of clogging of filters due to residual unsettled suspended matter.

Other measurements relate to the pH of water after chemical treatment (lime raises the value while iron and aluminium salts lower it) and the temperature of water at which the experiments are performed.

References

American Society for Testing Materials Norm ASTM D-2035-G4 T

Passino R. and M.Beccari (1970).

Standardization of jar test in coagulation-flocculation processes for the removal of inorganic turbidities (in Italian). Acqua e Aria, Milano, 1-10, October 1970.

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U.S. Association of Environmental Engineering Professors (1972). Environmental Engineering Unit Operations and Unit Processes Laboratory Manual. J.T. O'Connor Ed. 350 pp.

Leaching test

Leaching tests are used to simulate the behaviour of a waste in a dump that is submitted to the leaching activity of rain water. There are two different methods that use leaching solutions of different type:

- a) acetic acid, proposed for dumps receiving both organic and inorganic wastes;
- b) carbon dioxide saturated water, for dumps receiving only inorganic wastes.

Waste sample to be examined

If a liquid phase is present it must be removed by filtration or centrifugation and preserved in a refrigerator at 2-4°C. The granulometry of the solid phase, must be lower than 9.5 mm (standard sieve).

Waste extraction

20-25 g of the solid phase are weighed with a precision of ± 0.1 g.

Method a. A volume of distilled water corresponding to 16 times the weight of the sample is added.

The pH of the mixture is kept to 5 ± 0.2 using 0.5 N acetic acid. The stirring is performed for 24 hours, controlling the pH every 15 minutes during the first hour and then every hour.

The pH is corrected to 5 ± 2 using again 0.5 N acetic acid. If after 24 hours of stirring the pH is higher than 5.2, a new correction is made and the stirring prolonged for other 4 hours.

Method b. The leaching solution is given by a volume corresponding to 20 times the sample weight of distilled water saturated by carbon dioxide for not less than 15 minutes. The stirring is performed during 6 hours without any pH correction.

Stirring speed

Must be able to maintain continuously mixed the solid phase with the leaching solution.

Analysis of extract

The leaching solution and the liquid phase if present in the original sample are filtered by a membrane filter of 0.45 μm porosity and submitted to analysis for the single metals by the methods commonly used for waters.

GB Calculation of velocity gradient (G)

Formula

Velocity gradient (s^{-1})

$$G = \sqrt{P/\mu V}$$

Power requirement (W)

$$P = k * \rho * n^3 * D^5 \text{ turbulent}$$

$$P = k * \mu * n^2 * D^3 \text{ laminar}$$

Notes


Dynamic viscosity:	μ	0,001002	($N*s*m^{-2}$)
Mass density of fluid:	ρ	998,2	($kg*m^{-3}$)
Water volume:	V	0,001	(m^3)
Liquid mixed	Water		
Measurement temperature:	20		($^{\circ}C$)
Constant	k	3,38	
Number of turns / minute (RPM)	n	from 30 to 200	(min^{-1})
Number of turns / second	n	(from 30 to 200)*1/60	(s^{-1})
Diameter of impeller	D	0,075	(m)

Example: Velocity gradient for turbulent flow

Data

Number of turns / min-1	G (turbulent) / s^{-1}
200	544,28
150	353,52
100	192,43
75	124,99
50	68,03
40	48,68
30	31,62

Graph

GB Declaration of conformity

The manufacturer AQUALYTIC®
Adress Schleefstraße 12
 44287 Dortmund
 Germany

under our responsibility declare that the product is manufactured in conformity with the following standards:

EN 61010-1 (2001)

EN 61326-1 (1997) + A1 (1998) + A2 (2001) + A3 (2003)

2002/95/EC (RoHS)

2002/96/EC (WEEE)


and satisfies the essential requirements of the following directives:

Machines directive 2006/42/EC

Low voltage directive 2006/95/EC

Electromagnetic compatibility directive 2004/108/EC

plus modifications and that the documents listed in annex V are available at Velp's offices as foreseen by the machine directive.

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