## INSTRUCTION MANUAL: BOD TEST SET (Cat.-No. 418328)

### *Principle of Operation*

The calibration test can be used to check the operability of the whole OxiDirect BOD measurement system. A tablet is used for absorption of a defined amount of oxygen in the test system. This creates a reproducible negative pressure which can be used for a system check.

### *The operation is similar to a normal BOD measurement:*

- Set the temperature of the thermostatically controlled cabinet to  $20\text{ °C} \pm 0,1\text{ °C}$  and switch on, place the stirring system in the cabinet.
- Clean the bottles & stirring rods, use the overflow measuring flask to fill **157 ml** of distilled water in each bottle, add the stirring rod.
- Place the bottles in the thermostatically controlled cabinet without the sensor heads;
- After pre-tempering the distilled water filled sample bottles to  $20\text{ °C}$  place one test tablet in each bottle.
- Following please screw the sensor heads onto the bottles. The quiver is used as a gasket.
- Do not place KOH in the quiver!
- Then the measurement can be started (Measuring Range 0-400 mg/l; sample size 157ml)
- The Auto-start option should be switched off (ref. chapter 4.10.2)

### *Interpretation*

Compare the measured values against the check value. The measured value at the end of the 5-day test period should be within the tolerance range of the check value (the check value & tolerance are shown on the label on the packaging). Any differences between batches have already been taken into consideration. The check value will be obtained only if instructions are followed precisely. The highest test value is almost obtained after 24 hours. The negative pressure in the measuring system should then remain constant within the tolerance range till the end of the test period.

If only slight differences are noted between the measured value and the tolerance on the check value, the test should be repeated under precisely controlled conditions.

If significant differences occur between the measured value and the tolerance on the check value, the source of error may be determined as follows:

If the measured value decreases in the course of the test period, there is probably a leak in the system (damaged bottle thread, worn-out case). Replace the defective part and repeat the test. If the measured value remains constant outside the tolerance on the check value, several causes are possible: incompletely transferred tablet, contamination, wrong-sized bottle and/or test volume, wrong temperature in incubator, operating error or defective sensor.

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