
CHECKIT® micro pH-WP1 CHECKIT® micro pH-WP2

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
CHECKIT® micro pH-WP1 CHECKIT® micro pH-WP2

Conditioning

Remove the cap from the bottom of the tester to expose the electrode bulb and reference junction. Do not be alarmed if white crystals are present on the cap or electrode assembly. This is normal and it will dissolve during the conditioning pre-soak. For the conditioning pre-soak, use a small amount of tap water in a small cup and soak the electrode for at least 30 minutes. If the electrode dries out between uses, re-condition by using the above soaking procedure.

Calibration

Calibration is necessary and should be done regularly. The frequency is application dependent and can be determined by trial and error.

1. Turn the CHECKIT® micro WP1/2 on by pressing the ON/OFF button.
2. Dip about 2 – 3 cm of the electrode into a pH 7 buffer for general purpose use, or into pH 4 buffer (for usage in acidic solutions), or into pH 10 buffer (for usage in alkaline solutions). Use pH 7 buffer for calibration if the tests will be in closer to pH 7 than pH 4 or pH 10.
3. Press the CAL button to enter the calibration mode. A "CA" will flash on the display and then a pH value close to the pH buffer value will flash repeatedly.
4. When the flashing value stabilizes (normally after 6 seconds), press the HOLD/CON button to confirm the calibration.
5. The display will show "CO" and switch back to a pH reading of the buffer.
6. Rinse the electrode with de-ionized or tap water after use.

Possible problems with calibration

1. Failing to press the HOLD/CON button to enter the calibration. Pressing the CAL button instead will stop the flashing and resume measuring mode, but will not enter the calibration and the meter will not show the buffer value or measure accurately.
2. Failure to allow the CHECKIT® micro pH-WP1/2 time to sample the pH buffer till stabilization prior to pressing the HOLD/CON button confirming the calibration and small errors may be introduced.
3. Failure to re-hydrate the dry electrode after it has dried out before attempting a calibration will give fluctuating readings while it re-hydrates in a buffer and this causes errors.

pH-Testing

1. Remove cap from the electrode assembly and press the ON/OFF button to turn the CHECKIT® micro pH-WP1/2 on.
2. Dip the electrode 2 – 3 cm into a container filled with 20 – 30 ml of the test solution.
3. Stir gently for 10 – 20 seconds and let the reading stabilize.
4. Note the pH reading or press the HOLD/CON button to freeze the reading for logging.
(To release held reading, press HOLD/CON again.)
5. Press ON/OFF button to turn off the tester. If you do not press a button for 8.5 minutes, the tester will automatically shut-off to conserve battery power and prolong the battery life.

Error messages

E1 Batteries are low and need replacement.

E2 Wrong or bad buffer value (out of range) or the electrode is failing.
Check buffers or replace electrode.

OR Over range signal from voltage on a solution, the electrode is not contracting the solution, or the electrode is failing.
Check solution or replace electrode.

Maintenance

1. Rinse the electrode with de-ionized or tap water after each measurement to extend its useful life.
2. Periodic soaks in warm pH 4 buffer will help remove any contaminants that ruin electrodes.
3. If possible, keep a small piece of cloth or sponge in the electrode cap, moisten with clean water or electrode storage solution (not de-ionized water), and close the cap over the electrode.
4. The useful life of a CHECKIT® micro pH-WP1/2 is entirely dependent on the care the electrode and meter get.

The testers WP1 and WP2 are designed for water testing. Using them in other liquids may cause damage to the electrodes.

Changing batteries

1. Under the cap with the lanyard loop is the battery compartment. Wrap the rubber extrusion band provided around the battery compartment.
2. With dry hands firmly grip the cap and tester body. Twist the cap in anti-clockwise direction to open it.
3. Once it is opened, remove old batteries and replace with new ones noting the polarity. Replace the cap and twist tightly in clockwise direction until the black O-ring is under the cap and lanyard loop is lined with the face plate.
4. The cap must be tightened to ensure water tightness and it will require effort to tighten or loosen.

Electrode replacement

One of the benefits for selecting the waterproof CHECKIT® micro pH-WP1/2 now replace the meter portion of the instrument, thus saving money. When the tester fails to calibrate, gives fluctuating readings in buffers, or shows error messages **E2** or **OR** in a buffer solution, and the procedures in the maintenance section does not help, it is time to change the electrode.

Using the band provided, wrap it around the ribbed collar of the electrode module and slowly twist anti-clockwise to loosen. Twist until the collar is completely free and pull the electrode module straight out from the end of the tester. Keep and save both the ribbed collar and O-rings in the collar.

Align the new module so that its small slot on the CHECKIT® micro pH-WP1 and WP2 first (noting the connector pins on the electrode module) and gently slide the large tab into the large slot on the tester.

Gently push the electrode module onto the end to fully seat it in position.

Position the small O-ring inside the ribbed collar, and push the collar straight until it grips onto the first thread.

Twist is firmly in a clockwise direction until the larger O-ring that is close to the face plate disappears under the ribbed electrode module collar.


Important note

The CHECKIT® micro pH-WP1 and WP2 are designed for water testing. Using them in other liquids may cause damage to the electrode (see chapter "Maintenance").


Specifications WP1 / WP2

Model	WP1	WP2
pH-Range	0 to 14 pH	0 to 14 pH
Resolution	0.1 pH	0.1 pH
Accuracy	± 0.2 pH	± 0.1 pH
Calibration	4, 7 or 10	4, 7 und 10
Operating Temp.	0 to 50 °C	32 to 122 °F
ATC (Automatic Temperature Compensation)	no	yes
Battery Life	3 x 1.4 V	3 x 1.4 V
	Typical life = 60 hrs.	
Size (tester only)	ca. 165 x 38 mm	~
Weight (tester only)	ca. 90 gms	ca. 90 gms

Holding the CHECKIT® micro pH-WP

Schematic of sensor module

Re-order codes

1 Replaceable electrode for CHECKIT® micro pH-WP1	19 46 35
1 Replaceable electrode for CHECKIT® micro pH-WP2	19 46 45
100 Buffertablets pH 7	51 56 10
100 Buffertablets pH 4	51 56 20
100 Buffertablets pH 10	51 56 00
250 ml Deionised water	45 70 22
100 ml Plastic beaker	38 48 01
3 Replacement battery	33 32 09

pH-Buffer

Buffer solutions produced with tablets are colour coded for easy identification. The colour coding is as follow:

pH 4 tablet	= Red solution
pH 7 tablet	= Yellow solution
pH 10 tablet	= Green solution

Procedure

1. Fill a clean beaker with 20 ml deionised/distilled water.
2. Add one buffer tablet of the required value and shake or crush with a clean stirring rod until the tablet disintegrates.
3. Allow to stand for about two minutes to ensure that the buffer is completely dissolved and to allow any insoluble particles to settle.
4. Check that the colour of the solution corresponds to the buffer value required.
5. Immerse the electrode into the tube and calibrate the CHECKIT® micro in accordance with this instruction manual.
6. Remove the electrode and rinse it thoroughly in deionised/distilled water.

Note

It is important to rinse the electrode before immerse in a buffer solution of a different value.

Accuracy

Fresh buffer solutions should always be made up in deionised/distilled water.

The maximum tolerance based on the pH values 4, 7 and 10 is ±0.05 pH

Buffer Values

The pH value of buffer solutions varies with the temperature of solution.

The precise buffer values given by the buffer tablets at different solution temperatures are shown in the following table:

Temperature	Buffer value / pH		
	°C	°F	
0	32	4.01	7.12 10.33
5	41	4.01	7.09 10.25
10	50	4.00	7.06 10.18
15	59	4.00	7.04 10.11
20	68	4.00	7.02 10.05
25	77	4.01	7.00 10.00
30	86	4.01	6.99 9.95
35	95	4.02	6.98 9.92
40	104	4.03	6.98 9.88
45	113	4.04	6.97 9.85
50	122	4.06	6.97 9.82
55	131	4.07	6.98 9.80
60	140	4.09	6.98 9.77
70	158	4.12	6.99 9.73
80	176	4.16	7.00 9.69
90	194	4.20	7.02 9.66
Solution Colour	Red	Yellow	Green

Warranty

All Lovibond® CHECKIT® micros are guaranteed against defects in material or workmanship for a period of 6 months from the date of purchase. Please note however that the instrument is not guaranteed against misuse.

We reserve the right to alter or amend any of the items contained herein without prior notice.

Printed in Germany 10/96

Lovibond®, Tintometer® and CHECKIT® are registered trademarks of the Tintometer Group of Companies.


The Tintometer Ltd
Waterloo Road
Salisbury, SP1 2JY
Tel.: (+44) 1722 / 32 72 42
Fax: (+44) 1722 / 41 23 22
Telex: 4 73 72

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