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
Operation instructions

Content

1. General description	13
2. Principle of operation	13
3. Initial setup	14
3.1. Preperation for use	14
3.2. Electrical connection	14
4. Front-panel/operation controls and indicators	15
5. Operation	16
6. Service and maintenance	17
7. Troubleshooting	18
8. Technical data	19
9. Wiring diagram	20

1. General description

The Laboratory-Floc-Testers, models ET 710 and ET 720 with 4 and 6 stirring places are used as floc testers or universal agitators in laboratory and industrial environments.

The stirring rate is infinitely variable between 20 - 300 rpm. The stirring rate is shown in the LED display.

The stirring time is also set and shown via an LED display.

The high position of the agitator blades can be adjusted during operation, the rear panel is illuminated to allow a perfect observation of the samples.

2. Principle of operation

The instrument is driven by one common motor for all stirring positions. This ensures a uniform stirring rate for all stirring positions.

The individual stirring positions are connected with the motor by means of pulley belts which can be replaced after removal of the upper covering plate. The speed of the motor is directly calculated with the aid of an incremental transducer and after electronic evaluation transmitted to the revolution indicator.

Presetting of the stirring time is made possible through the timer-LED-display. The stirring time can be entered in increments of one minute each. When this preset time has elapsed the stirrer is switched-off automatically.

The shafts of the agitator blades are fitted with slip clutches close to the bearings, this allows the manual stopping of individual stirring positions as well as high adjustments during operation.

Attention! Regardless of the slip clutches it is dangerous to reach into the stirring position whilst the mixing blades are still revolving. Risk of injury! Do not operate the instrument unattended. Operation by qualified personnel only and after careful review of this instruction manual.

3. Initial setup

Unpacking the instrument

When removing the instrument from the shipping box inspect instrument and accessories for any damage that may have occurred during the shipment. To assure entitlement to insurance cover, all damages must be reported to the carrier immediately.

Verification of the shipment

Verify that the following items are present:

- models ET 710 / ET 720 (4 or 6 places)
- power cord
- instrument manual
- rubbermat (for the beakers)

Storage

The instrument requires no special treatment, it should, however, be packed impactproof.

Avoid high humidity during storage!

3.1. Preperation for use

Place the instrument on a laboratory bench for easy accessibility and operation. In the immediate neighbourhood of the instrument no high current consuming apparatus (such as motors etc.) may be operated.


3.2. Electrical connection

Attention! Opening of the instrument may be performed by qualified personnel only. The VDE guidelines must be observed.

Compare the available line voltage with the voltage required for the instrument (see name plate). Defects resulting from the non-observation of this instruction are not covered by the warranty. The electric power supply must be provided with an operating ground-connection.

If line voltage and voltage rating of the instrument are identical set power switch (1) to "0" and connect the power cord.

4. Front-panel/operation controls and indicators

- ① Power switch
- ② Switch for rear panel illumination
- ③ Display stirring rate range
- ④ Increase stirring rate
- ⑤ Decrease stirring rate
- ⑥ Display timer
- ⑦ Timer switch-on
- ⑧ Timer switch-off
- ⑨ Increase timer
- ⑩ Decrease timer

5. Operation

- Place the beakers containing the samples under the respective agitator blades and lower the agitators to the required position. To guard against inadvertent contact with the agitator blades always fit-out each available stirring position with a beaker (eventually without sample) even if it is not required for the test.
- Switch-on the instrument by setting the power switch (1) (see figure on page 5) in position "I".
- Set the switch for the back panel-illumination (2) in position (1), if illumination of the sample is desired.

Changing the stirring rate-setting

- By means of the push buttons (4) and (5) the rate can be set in the range between 20 and 300 revolutions per minute. To do this keep the button pressed down until the desired rate appears in the display. After setting wait for the hunting until the rate display has stabilized. A "000" display means that the motor does not rotate.

Attention! If the "000" is underpassed when setting the stirring rate the control (ring-mode) switches to 300 rpm.

Changing the timer-setting

- By pressing push button (7) the display of the timer is switched on. The standard setting of 1 minute is displayed.
- By means of the push buttons (9) and (10) the time can be set in the range between 1 minute and 999 minutes. The minimum increment is 1 minute.
- Pressing push button (7) again activates the timer and flashes until the set time has elapsed.
- When the set time has elapsed the flashing stops, the time remains however as set and the rate display goes to "000". Pressing push button (7) again initiates another stirring sequence of the same duration.
- By means of push button (7) the timer can also be switched-off during operation or it can be switched-over to continuous operation by means of push button (8).

6. Service and maintenance

The instrument needs no special maintenance or care. If the front panel of the instrument is dirty, it can be cleaned with a moist tissue. **Be careful that no water enters the instrument !**

After each test the agitator blades should be cleaned with water to avoid any risk of contamination or modification of the next sample.

Stopsetting / not in use

In case of stopsetting or if the instrument is not of be used for a longer period generally no special precautions are required.

Storage

If the instrument is to be stored, no special precautions are required.

For storage purposes, however, shock proof packing and protection from moisture is recommended.

Preparation for transport

The instrument must be dry for the transport. Dry it beforehand if necessary.

The instrument must be transported in shock-proof packing (original packing, if possible). The packing must be marked accordingly (keep dry/fragile).

7. Troubleshooting

Attention! Before opening the instrument disconnect the power cord from the mains outlet. Never work on a voltage-carrying instrument. Remount the instrument completely after any repair prior to testing it.

The following list should enable you to detect any malfunction immediately and to restore the correct function at reasonable cost. Please read the information about causes and possible consequences carefully and review the respective paragraphs of this manual for assistance.

One stirring place doesn't work, all other places function correctly !

- The pulley belt of this stirring place is broken. Before starting the replacement of this belt disconnect the power. Then unscrew the 4 or 6 handles of the stirring places. Remove the screws of the upper covering plate and take-off the plate. Replace the broken belt and remount the instrumente in reverse order.

None of the stirring places works !

- Pulley belt of the motor broken. Before starting the replacement of this belt disconnect the power. Repair as described above.

All stirring places function correctly, however no display of the rate !

- Incremental transducer on the motor is defect.
Replacement of this unit at the factory of the manufacturer only.
Please observe the instructions for the return of the instrument.

Both displays don't function after the instrument has been switched-on !

- Check the fuse on the back plate under the instrument jack (power cord receptacle) slow-blow fuse 1 A/250 V (2 x).

Both displays don't function after the instrument has been switched-on !

- The circuit board of the display or of the control is defect. Replacement of either unit at the factory of the manufacturer only. Please observe the instructions for the return of the instrument.

Instruction for the return of the instrument

If repair of the a.m. malfunction should prove to be impossible (review the operating instructions to exclude any improper operation) the instrument must be returned. The repair time for instruments sent back without an enclosed and detailed description of the malfunction may be delayed, some defects might remain undiscovered which in turn can lead to an incomplete repair.

The instrument must be packed shock-proof - preferably in the original shipping box. See also paragraph **6 Service and maintenance**

8. Technical data

Power requirement:	230 V / 50 A
Power consumption:	150 VA
Fuse:	2 x 1 A slow-blow, 250 V, 5 x 20 mm
Protective system:	IP 31
Enclosure:	painted iron steel
Stirring rate:	20 - 300 revolutions / minute
Accuracy:	± 10 revolutions / minute
Timer:	1 - 999 minutes
Resolution:	1 minute
Ambient temperature:	10 °C - 45 °C


ET 710

weight:	23 kg
Dimensions stirring compartment:	250 h x 650 w x 150 d mm
Dimensions enclosure:	460 h x 700 w x 240 d mm
Art.No.:	2 41 91 30

ET 720

weight:	31 kg
Dimensions stirring compartment:	250 h x 950 w x 150 d mm
Dimensions enclosure:	460 h x 990 w x 240 d mm
Art.No.:	2 41 91 31

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We reserve the right to alter or amend
any of the items contained herein without
prior notice.

Lovibond® is a registered trademarks of
the Tintometer Group of Companies.

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We reserve the right to alter or amend
any of the items contained herein
without prior notice.

AQUALYTIC GmbH & Co.
Postfach 1143
D-63231 Neu-Isenburg
Germany

Telefon +49 (0) 61 02 - 78 94 - 5
Telefax +49 (0) 61 02 - 78 94 99