MiniDirect Photometer

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CLt

Operation

Switch the unit on using the ON/OFF key.

The display shows the following:

Attention:

The unit is designed to measure chlorine either with tablet reagents (CLt) or with powder packs (CLP). Please select the required chlorine method as described below.

Mode

METHOD

Select the test required using the MODE key:

CLt → CLP → PH → .CYS → CLt → (Scroll)

The display shows the following:

Fill a clean vial with the water sample up to the 10 ml mark, screw the cap on and place in the sample chamber with the Δ-mark on the vial aligned with the ∇-mark on the instrument.

Zero
Test

METHOD

0.0.0

Press the ZERO/TEST key.

The method symbol flashes for approx. 8 seconds.

The display shows the following:

After zero calibration is completed, remove the vial from the sample chamber. Add the appropriate reagent; a colour will develop in the sample.

Screw the cap back on and place the vial in the sample chamber with the Δ and ∇ marks aligned.

Zero
Test

METHOD

RESULT

Press the ZERO/TEST key.

The method symbol flashes for approx. 3 seconds.

The result appears in the display.

The result is saved automatically.

Repeating the analysis:

Press the ZERO/TEST key again.

New zero calibration:

Press the MODE key until the desired method symbol appears in the display again.

Zero
Test

Mode

Display backlight

Press the "!" key to turn the display backlight on or off. The backlight is switched off automatically during the measurement.

!

Technical Data

Light source:	LED, Filter ($\lambda = 528 \text{ nm}$)
Battery:	9 V-block battery (Life 600 tests), without display light.
Auto-OFF:	Automatic switch off 5 minutes after last keypress
Ambient conditions:	5-40°C 30-90 % rel. humidity (non-condensing).
CE:	Certificate for Declaration of CE-Conformity at www.aqualytic.de


Method notes

Prior to measurement ensure that the sample is suitable for analysis (no major interferences) and does not require any preparation i.e. pH adjustment, filtration etc.


Reagents are designed for use in chemical analysis only and should be kept well out of the reach of children. Ensure proper disposal of reagent solutions.

Material Safety Data Sheets: www.aqualytic.de

Correct filling of the vial

Opening powder packs

Chlorine 0.01 - 6.0 mg/l with tablets (CLt)

CLt

0.0.0

Zero
Test

CLt

RESULT

(a) free chlorine

Perform zero calibration (see "Operation").

Empty the vial and then add a DPD No. 1 tablet. Crush the tablet with a clean stirring rod then add the water sample to the 10 ml mark. Mix well with the stirring rod to dissolve the tablet. Screw the cap on and replace the vial in the sample chamber making sure the Δ and ∇ marks are aligned.

Press the ZERO/TEST key.

The method symbol flashes for approx. 3 seconds.

The result is shown in the display in mg/l free chlorine.

(b) total chlorine

Remove the vial and add one DPD No. 3 tablet to the coloured test solution. Mix to dissolve with the stirring rod. Replace the cap and put the vial back into the sample chamber, repositioning the Δ and ∇ marks.

Wait for a colour reaction time of two minutes.

Press the ZERO/TEST key.

The method symbol flashes for approx. 3 seconds.

The result is shown in the display in mg/l total chlorine.

Rinse the vial and cap thoroughly after each test.

Zero
Test

CLt

RESULT

(c) combined chlorine

combined chlorine = total chlorine - free chlorine

Measuring tolerance: 0-1 mg/l: ± 0.05 mg/l > 3-4 mg/l: ± 0.30 mg/l
 > 1-2 mg/l: ± 0.10 mg/l > 4-6 mg/l: ± 0.40 mg/l
 > 2-3 mg/l: ± 0.20 mg/l

Chlorine 0.01 - 2.0 mg/l with Powder Packs (CLP)

CLP

0.0.0

Zero
TestCLP

RESULT

(a) free chlorine

Perform zero calibration (see "Operation").

Remove the vial from the sample chamber. Add a VARIO Chlorine FREE-DPD/F10-powder pack. Screw the cap on and swirl the vial for 20 seconds. Replace the vial in the sample chamber immediately making sure the Δ and ∇ marks are aligned.

Press the ZERO/TEST key.

The method symbol flashes for approx. 3 seconds.

The result is shown in the display in mg/l free chlorine.

Rinse the vial and cap thoroughly after each test.

(b) total chlorine

Perform zero calibration (see "Operation").

Remove the vial from the sample chamber. Add a VARIO Chlorine TOTAL-DPD/F10-powder pack. Screw the cap on and swirl the vial for 20 seconds. Replace the vial in the sample chamber making sure the Δ and ∇ marks are aligned.

Wait for a colour reaction time of 3 minutes.

Press the ZERO/TEST key.

The method symbol flashes for approx. 3 seconds.

The result is shown in the display in mg/l total chlorine.

Rinse the vial and cap thoroughly after each test.

Zero
TestCLP

RESULT

(c) combined chlorine

combined chlorine = total chlorine - free chlorine

Measuring tolerance: 0-1 mg/l: ± 0.05 mg/l

> 1-2 mg/l: ± 0.10 mg/l

pH value 6.5 - 8.4 with tablets (PH)

PH

0.0.0

Zero
Test

PH

RESULT

Perform zero calibration (see "Operation").

Remove the vial from the sample chamber. Add a PHENOLRED/PHOTOMETER tablet and mix to dissolve using a clean stirring rod. Screw the cap on and replace the vial in the sample chamber making sure the Δ and ∇ marks are aligned.

Press the ZERO/TEST key.

The method symbol flashes for approx. 3 seconds.

The pH value is shown in the display.

Rinse the vial and cap thoroughly after each test.

Measuring tolerance: ± 0.1 pH

pH value 6.5 - 8.4 with liquid reagent (PH)

PH

0.0.0

Zero
Test

PH

RESULT

Perform zero calibration (see "Operation").

Hold the reagent bottle vertically and add 6 evenly sized drops of PHENOLRED-solution to the vial by pressing gently. Replace the cap tightly, swirl and place the vial in the sample chamber making sure the Δ and ∇ marks are aligned.

Press the ZERO/TEST key.

The method symbol flashes for approx. 3 seconds.

The pH value is shown in the display.

Rinse the vial and cap thoroughly after each test.

Measuring tolerance: ± 0.2 pH

GB

Mode

On
Off

!

Mode

Menu selections

Press the MODE key and hold.

Switch the unit on using the ON/OFF key.
Allow the 3 decimal points to be displayed before releasing the MODE key.

The "!" key allows for selection of the following menu points:

- recall stored data
- setting the date and time
- user calibration

The selected menu is indicated by an arrow in the display.

Confirm the selection with the MODE key.

Recall of stored data

The meter shows the most recent measurements taken in the following format (automatically proceeds every 3s until result is displayed):

Number n xx (xx: 16...1)
 Year YYYY (i.e. 2007)
 Date mm.dd (monthmonth:dayday)
 Time hh:mm (hourhour:minute) minute
 Test Method
 Result x,xx

Zero
Test

Mode

!

The ZERO/TEST key repeats the current data set.

The MODE key scrolls through all stored data sets.


Quit the menu by pressing "!" key.

Setting date and time (24-hour-format)

After confirming the selection with the MODE key the value to be edited will be shown for 2 sec.

The setting starts with the year (YYYY) followed by the actual value to be edited. Same applies for month (mm), day (dd), hour (hh) and minutes (mm). Set the minutes first in steps of 10, press the "!" key to continue setting of minutes in steps of 1.

Mode

Mode

Zero
Test

!

Increase the value by pressing the MODE key.

Decrease the value by pressing ZERO/TEST key.

Proceed to the next value to be edited by pressing "!" key.

After setting the minutes and pressing the "!" key the display will show "IS SET" and instrument returns into the measurement mode.

When the battery is taken of for more than 1 minute, the unit will automatically enter the date/time menu when switched on again.

User calibration

Note:

cAL

user calibration (Display in calibration mode)

CAL

factory calibration (Display in calibration mode)

Mode

On
Off

!

CAL

CLt

Press MODE key and **hold**.

Switch the unit on using the ON/OFF switch,
Release the MODE key after approx. 1 second.
Select the required menu points using the "!" key.

After confirming the selection with the MODE key the instrument will show CAL/CLt.

Scroll through methods using the MODE key.

Perform zero calibration (see "Operation").

Press the ZERO/TEST key.

The method symbol flashes for approx. 8 seconds.

The display shows the following in alternating mode:

Perform calibration with a standard which concentration is known
(see "Operation").

Press the ZERO/TEST key.

The method symbol flashes for approx. 3 seconds.

The result is shown in the display, alternating with CAL.

If the reading corresponds with the value of the calibration standard (within the specified tolerance), exit calibration mode by pressing the ON/OFF key.

Changing the shown value:

Otherwise, pressing the MODE key once increases the displayed value by 1 digit.

Pressing the ZERO/TEST key once decreases the displayed value by 1 digit.

Press the corresponding key until the reading equals the value of the calibration standard.

By pressing the ON/OFF key, the new correction factor is calculated and stored in the user calibration software.

Confirmation of calibration (3 seconds).

Zero
Test

METHOD

0.0.0

CAL

Zero
Test

METHOD

RESULT

CAL

Mode

Zero
Test

CAL

RESULT + x

On
Off

: :

Recommended calibration values

Chlorine: between 0.5 and 1.5 mg/l

pH: between 7.6 and 8.0

Cyanuric acid: between 30 and 60 mg/l

Resetting the factory calibration

Resetting the user calibration to the original factory setting resets all parameters.

A user calibrated method will be indicated by an arrow in the display.

To reset the calibration to the factory setting:

Press both the MODE and ZERO/TEST key and **hold**.

Switch the unit on using the ON/OFF key. Release the MODE and ZERO/TEST keys after approx. 1 second.

The following messages will appear in turn on the display:


The factory setting is active.
(SEL stands for Select)

or:


Calibration has been set by the user.
(If the user calibration is to be retained, switch the unit off using the ON/OFF key.)

Calibration is reset to the factory setting by pressing the MODE key.
The following messages will appear in turn on the display:

Switch the unit off using the ON/OFF key.

Changing the battery

Error codes

General error codes

EOOI	Light absorption too great. Reasons: dirty optics.
HI	Measuring range exceeded or excessive turbidity.
LO	Result below the lowest limit of the measuring range.
LO BAT	Replace 9 V battery, no further tests possible.
EOIO	Calibration factor "out of range"

Individual error codes

E 70	CLt: Factory calibration incorrect / erase
E 72	CLP: Factory calibration incorrect / erase
E 74	PH: Factory calibration incorrect / erase
E 76	.Cys: Factory calibration incorrect / erase
E 71	CLt: User calibration incorrect / erase
E 73	CLP: User calibration incorrect / erase
E 75	PH: User calibration incorrect / erase
E 77	.Cys: User calibration incorrect / erase

Tintometer GmbH, Division Aqualytic®

Schleefstraße 8-12
D-44287 Dortmund
Tel.: (+49) (0)2 31 / 9 45 10-755
Fax: (+49) (0)2 31 / 9 45 10-750
sales@aqualytic.de
www.aqualytic.de
Germany

Technical changes without notice
Printed in Germany 11/09