MiniDirect Photometer

Eisen • Iron

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(GB) Instruction Manual

DE

On
Off

Fet

Mode

METHODE

Zero
Test

METHODE

0.0.0

Zero
Test

METHODE

ERGEBNIS

Zero
Test

Mode

!

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Operation

Switch the unit on using the ON/OFF key.

The display shows the following:

Attention:

The unit is designed to measure iron either with tablet reagents (Fet) or with powder packs (Fe:P1/Fe:P2). Please select the required chlorine method as described below.

Select the test required using the MODE key:

Fet → Fe:P1 → Fe:P2 → Fet → (Scroll)

The display shows the following:

Fill a clean vial with the water sample up to the 10 ml mark, screw the cap on and place in the sample chamber with the Δ-mark on the vial aligned with the ∇-mark on the instrument.

Press the ZERO/TEST key.

The method symbol flashes for approx. 8 seconds.

The display shows the following:

After zero calibration is completed, remove the vial from the sample chamber. Add the appropriate reagent; a color will develop in the sample.

Screw the cap back on and place the vial in the sample chamber with the Δ and ∇ marks aligned.

Press the ZERO/TEST key.

The method symbol flashes for approx. 3 seconds.

The result appears in the display.

The result is saved automatically.

Repeating the analysis:

Press the ZERO/TEST key again.

New zero calibration:

Press the MODE key until the desired method symbol appears in the display again.

Display backlight

Press the "!" key to turn the display backlight on or off. The backlight is switched off automatically during the measurement.

Technical Data

Light source:	LED, Filter ($\lambda = 528 \text{ nm}$)
Battery:	9 V-block battery (Life 600 tests), without display light.
Auto-OFF:	Automatic switch off 10 minutes after last keypress
Ambient conditions:	5-40°C 30-90 % rel. humidity (non-condensing).
CE:	Certificate for Declaration of CE-Conformity at www.aqualytic.de


Method notes

Prior to measurement ensure that the sample is suitable for analysis (no major interferences) and does not require any preparation i.e. pH adjustment, filtration etc.


Reagents are designed for use in chemical analysis only and should be kept well out of the reach of children. Ensure proper disposal of reagent solutions.

Material Safety Data Sheets: www.aqualytic.de

Correct filling of the vial

Opening powder packs

GB

Iron (II- and III-Ions) 0.02 - 1.0 mg/l Fe with tablets (Fet)

0.0.0

Zero
Test

Fet

RESULT

Perform zero calibration (see "Operation").

Add one IRON LR tablet straight from the foil to the 10 ml water sample, and crush using a clean stirring rod. Mix well with the stirring rod to dissolve the tablet. Screw the cap on and replace the vial in the sample chamber making sure the Δ and ∇ marks are aligned.

Wait for a colour reaction time of 5 minutes!

Press the ZERO/TEST key.

The method symbol flashes for approx. 3 seconds.


The result is shown in the display in mg/l total dissolved iron.

Measuring tolerance: $\pm 0,05$ mg/l

Notes

The IRON (II) LR tablet is used for differentiation - as described above - in place of the IRON LR tablet.

Differentiation

* Supplementary Procedure

Add 1 ml of concentrated sulphuric acid to 100 ml of the water sample. Heat and boil for 10 minutes or until all particles have dissolved. After cooling down the sample is set on a pH-value of 3 to 6 by using ammonia solution. Refill with distilled water to the previous volume of 10 ml. Mix well. Pour into the vial and fill to the 10 ml mark. Add an IRON LR-tablet, crush and mix well to dissolve. Allow to stand for 5 minutes. Water which has been treated with organic compounds as corrosion inhibitors must be oxidised where necessary to break down the iron complexes - add 1 ml of concentrated sulphuric acid and 1 ml of concentrated nitric acid to a 100 ml sample and boil to approximately half volume. After cooling down proceed with the analysis as described above.

** not included

Iron 0.02 - 3 mg/l Fe with Powder Packs (Fe:P1)

0.0.0

Zero
Test

Fe:P1

RESULT

Perform zero calibration (see "Operation").

Add **one** Vario Ferro F10 powder pack straight from the foil to the water sample.

Close the vial with the cap tightly and swirl the vial several times to mix the contents (Note 4).

Place the vial in the sample chamber making sure that the Σ marks are aligned.

Wait for a **colour reaction time of 3 minutes**

Press ZERO/TEST key.

The method symbol flashes for approx. 3 seconds.

The result is shown in the display in mg/l iron.

Measuring tolerance: $\pm 0,05$ mg/l Fe.

Notes

1. The reagent reacts with all soluble iron and most insoluble forms of iron in the water sample.
2. Iron oxide requires a prior digestion, use mild, vigorous or Digesdahl digestion (e.g. for digestion with acid see page 91).
3. Very strong alkaline or acidic water samples must be adjusted to a pH-Value between 3 and 5 before analysis.
4. Accuracy is not affected by undissolved powder.
5. Water samples containing visible rust should be allowed to react at least five minutes.

Iron, total (TPTZ) 0.05 - 1.8 mg/l Fe with Powder Packs (Fe:P2)

0.0.0

Do **not** perform zero calibration.

Use two clean vials (24 mm Ø) and mark one as blank for zeroing.

Fill a clean vial with 10 ml deionized water (this is the blank).

Fill the second clean vial with 10 ml water sample (this is the sample).

Add **one** Vario IRON TPTZ F10 powder pack straight from the foil into each vial.

Close the vials with the caps tightly and swirl the vials several times to mix the contents.

Wait for a **colour reaction time of 3 minutes**.

After reaction period is finished proceed as follows:

Place the vial (the blank) in the sample chamber making sure that the Σ marks are aligned.

Press ZERO/TEST key.

The method symbol flashes for approx. 8 seconds.

Remove the vial from the sample chamber.

Place the vial (the sample) in the sample chamber making sure that the Σ marks are aligned.

Press ZERO/TEST key.

The method symbol flashes for approx. 3 seconds.

The result is shown in the display in mg/l Iron.

Measuring tolerance: $\pm 0,05$ mg/l Fe.

Zero
Test

Fe:P2

RESULT

Zero
Test

Fe:P2

RESULT

Notes

1. For determination of total Iron digestion is required.
TPTZ reagent recovers most insoluble iron oxides without digestion.
2. Rinse all glassware with 1:1 Hydrochloric acid solution first and then rinse with deionised water to remove iron deposits that can cause slightly high results.
3. Strong alkaline or acidic water samples must be adjusted to between pH 3 and pH 8 before the reagent is added (use 0.5 ml Sulfuric acid resp. 1 mol/l Sodium hydroxide).

GB

Mode

On
Off

!

Mode

Menu selections

Press the MODE key and hold.

Switch the unit on using the ON/OFF key.
Allow the 3 decimal points to be displayed before releasing the MODE key.

The "!" key allows for selection of the following menu points:

- ▲ recall stored data
- ▲ ▽ setting the date and time
- ▼ user calibration

The selected menu is indicated by an arrow in the display.

Confirm the selection with the MODE key.

Recall of stored data

The meter shows the most recent measurements taken in the following format (automatically proceeds every 3s until result is displayed):

Number n xx (xx: 16...1)
 Year YYYY (i.e. 2007)
 Date mm.dd (monthmonth:dayday)
 Time hh:mm (hourhour:minute) minute)
 Test Method
 Result x,xx

Zero
Test

Mode


!

The ZERO/TEST key repeats the current data set.

The MODE key scrolls through all stored data sets.

Quit the menu by pressing "!" key.

Mode

Mode

Zero
Test

!

Setting date and time (24-hour-format)

After confirming the selection with the MODE key the value to be edited will be shown for 2 sec.

The setting starts with the year (YYYY) followed by the actual value to be edited. Same applies for month (mm), day (dd), hour (hh) and minutes (mm). Set the minutes first in steps of 10, press the "!" key to continue setting of minutes in steps of 1.

Increase the value by pressing the MODE key.

Decrease the value by pressing ZERO/TEST key.

Proceed to the next value to be edited by pressing "!" key.

After setting the minutes and pressing the "!" key the display will show "IS SET" and instrument returns into the measurement mode.

When the battery is taken of for more than 1 minute, the unit will automatically enter the date/time menu when switched on again.

User calibration

Note:

cAL

user calibration (Display in calibration mode)

CAL

factory calibration (Display in calibration mode)

Mode

On
Off

!

CAL

Fet

Press MODE key and **hold**.

Switch the unit on using the ON/OFF switch,
Release the MODE key after approx. 1 second.
Select the required menu points using the "!" key.

After confirming the selection with the MODE key the instrument will show CAL/Fet.

Scroll through methods using the MODE key.

Perform zero calibration (see "Operation").

Zero
Test

Press the ZERO/TEST key.

METHOD

The method symbol flashes for approx. 8 seconds.

0.0.0

The display shows the following in alternating mode:

CAL

Perform calibration with a standard which concentration is known
(see "Operation").

Zero
Test

Press the ZERO/TEST key.

METHOD

The method symbol flashes for approx. 3 seconds.

RESULT

The result is shown in the display, alternating with CAL.

CAL

If the reading corresponds with the value of the calibration standard (within the specified tolerance), exit calibration mode by pressing the ON/OFF key.

Changing the shown value:

Mode

Otherwise, pressing the MODE key once increases the displayed value by 1 digit.

Zero
Test

Pressing the ZERO/TEST key once decreases the displayed value by 1 digit.

CAL

Press the corresponding key until the reading equals the value of the calibration standard.

RESULT + x

On
Off

By pressing the ON/OFF key, the new correction factor is calculated and stored in the user calibration software.

:

Confirmation of calibration (3 seconds).

Recommended calibration values

Iron (Fet/Fe:P2): between 0.3 and 0.7 mg/l Fe

Iron (Fe:P1): between 0.1 and 2.0 mg/l Fe

Resetting the factory calibration

Resetting the user calibration to the original factory setting resets all parameters.

A user calibrated method will be indicated by an arrow in the display.

To reset the calibration to the factory setting:

Press both the MODE and ZERO/TEST key and **hold**.

Switch the unit on using the ON/OFF key. Release the MODE and ZERO/TEST keys after approx. 1 second.

The following messages will appear in turn on the display:

SEL

CAL

The factory setting is active.
(SEL stands for Select)

or:

SEL

cAL

Calibration has been set by the user.
(If the user calibration is to be retained, switch the unit off using the ON/OFF key.)

Mode

SEL


CAL

On
Off

Calibration is reset to the factory setting by pressing the MODE key.
The following messages will appear in turn on the display:

Switch the unit off using the ON/OFF key.

Changing the battery

Error codes

General error codes

E00I

Light absorption too great. Reasons: dirty optics.

HI

Measuring range exceeded or excessive turbidity.

LO

Result below the lowest limit of the measuring range.

LO BAT

Replace 9 V battery, no further tests possible.

E0IO

Calibration factor "out of range"

Individual error codes

E 70

Fet: Factory calibration incorrect / erase

E 72

Fe:P1: Factory calibration incorrect / erase

E 74

Fe:P2: Factory calibration incorrect / erase

E 71

Fet: User calibration incorrect / erase

E 73

Fe:P1: User calibration incorrect / erase

E 75

Fe:P2: User calibration incorrect / erase

Tintometer GmbH, Division Aqualytic®

Schleefstraße 8-12
D-44287 Dortmund
Tel.: (+49) (0)2 31 / 9 45 10-755
Fax: (+49) (0)2 31 / 9 45 10-750
sales@aqualytic.de
www.aqualytic.de
Germany

Technical changes without notice
Printed in Germany 11/09